

1. Device Name	POLY GLOVES																			
2. Description	Low Density Polythylene (LDPE) Disposable Glove, Non Medical Grade, Two Piece Construction, Single Heat Sealed Seam, Embossed Finish , Ambidextrous, Non Sterile, Single Use.																			
3. Product Codes & Sizes	<table><tr><th>Product Code</th><th>Size</th><th>Inner Barcode</th><th colspan="2">Shipper Barcode</th></tr><tr><td>41053#</td><td>MEDIUM</td><td>9311059410537</td><td colspan="2">19311059410531</td></tr><tr><td>41054#</td><td>LARGE</td><td>9311059410544</td><td colspan="2">19311059410548</td></tr></table>					Product Code	Size	Inner Barcode	Shipper Barcode		41053#	MEDIUM	9311059410537	19311059410531		41054#	LARGE	9311059410544	19311059410548	
Product Code	Size	Inner Barcode	Shipper Barcode																	
41053#	MEDIUM	9311059410537	19311059410531																	
41054#	LARGE	9311059410544	19311059410548																	
4. Recommended Applications	- Food Processing - Food Service - Food Manufacturing - Cafes - Hair and Beauty - Childcare - Cleaning - Hygienic - Sandwich Bars - Hospitality																			
5. Materials	Low density polyethylene film grade (complies with FDA 177.1520 for food contact)																			
6. Production Process	Materials, compounding, extrusion, heating (155-177c), filtering, cooling, rolling extrusion film, QC film gauge, Glove stamping with heat seal die cut, stripping gloves, QC inspection and tests, packing, final QC inspection and tests, Dispatch shipment.																			
7. Standards & Approvals																				



11. Quality Control Inspection Procedures	a) Shipment is placed in quarantine area separately from other gloves in warehouse as one lot b) Shipment is then sorted by sizes if shipment consists of more than one size. c) QC inspector records lot number(s) indicted on inner boxes and outer cases, sizes, quantity of shipment, buyer, date, etc d) QC inspector selects at random samples according to ISO procedures for sampling, each lots of different sizes are handled separately e) Samples are sent for testing as per requirements f) Results of tests will be recorded and only those that meet acceptance criteria are released for dispatch. If shipment consists of more than one size, shipment will be released only when all lots of different sizes have passed. g) If a shipment is found to fail test requirements it will be rejected by QC inspector and not released for shipment.															
12. Packaging	500 pieces will be packed into a dispenser box. The number contained in each dispenser box determined by count. 4 dispenser boxes are packed into 3 layer inner carton, 5 inner cartons per export quality shipper. For a total of 20 dispenser boxes per shipper carton.															
13. Shipper Carton Specifications	<table><tr><th>SIZE</th><th>LENGTH (mm)</th><th>WIDTH (mm)</th><th>HEIGHT (mm)</th><th>GROSS WEIGHT (kg)</th></tr><tr><td>MEDIUM</td><td>700</td><td>290</td><td>265</td><td>13.50</td></tr><tr><td>LARGE</td><td>700</td><td>290</td><td>265</td><td>14.0</td></tr></table>	SIZE	LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)	GROSS WEIGHT (kg)	MEDIUM	700	290	265	13.50	LARGE	700	290	265	14.0
SIZE	LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)	GROSS WEIGHT (kg)												
MEDIUM	700	290	265	13.50												
LARGE	700	290	265	14.0												
14. Storage Conditions	The gloves will maintain their properties when stored in dry conditions out of direct sunlight at between 10°c to 30°c.															
15. Shelf Life	Not Specified															

Poly Gloves

500 medium

polyethylene disposable gloves

powder free
ambidextrous

